



PRIYAN INTERNATIONAL LAB AND TECHNOLOGY

(Reference Material Product as per International Standard ISO: 17034)

2nd and 3rd Floor, C-247, Sector-10, Noida, Gautam Buddha Nagar, Uttar Pradesh-201301

Tel No. 0120-3684527, Mob.: 8882764797

Mail: priyanintlabetech@gmail.com, Web: www.priyaninternationalabtech.com


CRM CERTIFICATE

Format No: PILT/QSP/055/00/FMT/02

Certificate No: PILT/CRM/M-009/23/019


METHYL COBALAMIN

STRUCTURE:



DESCRIPTION & IDENTIFICATION:

Batch No.: PILTRS/23/01/019	Ref./Product No.: M- 009
Unit Quantity: 1 gm	Chemical Formula: $\text{C}_{63}\text{H}_{91}\text{CoN}_{13}\text{O}_{14}\text{P}$
Molecular Weight: 1344.40 g/mol	Assigned Value (Purity): 88.1 % w/w or 0.881 mg per mg on as is basis
Date of Release: 20/02/2024	uCRM(%) = 0.2209%
Validity Date: 19/02/2027	Method: IP 2022
Date of Issue: 20/06/2024	Storage: Keep container tightly closed, protected from light and store between 2°C to 8°C temp.



PLATINUM PLUS
ISO 17034
Accredited
Reference
Material Producer





PRIYAN INTERNATIONAL LAB AND TECHNOLOGY

(Reference Material Product as per International Standard ISO: 17034)

2nd and 3rd Floor, C-247, Sector-10, Noida, Gautam Buddha Nagar, Uttar Pradesh-201301

Tel No. 0120-3684527, Mob.: 8882764797

Mail: priyanintlabetech@gmail.com, Web: www.priyaninternationalabtech.com

CRM CERTIFICATE

Format No: PILT/QSP/055/00/FMT/02

UNCERTAINTY:

The assigned uncertainty covers uncertainty contribution from characterization, in homogeneity, storage & transport stability etc. (wherever applicable), is the combined standard uncertainty, calculated using a coverage factor ($K=2$) which gives a level of confidence of approx. 95%. As per ISO 17034:2016 & ISO Guide 35, for this pharmaceutical standard assigned uncertainty value is considered to be negligible w.r.t. defined limits of method specific assays for which the PILTRM/CRM is used.

METROLOGICAL TRACEABILITY AND MEASUREMENT METHODS:

NIST or other traceable standards are used for calibration and performance verification of instruments. The assigned value is traceable to SI units through the use of Primary Standard Mass Balance Methods (Physical and chemical). Characterization was done by the combination of Primary Reference Methods viz. NMR, LCMS, FTIR with use of pure substance/traceable RM/CRM in compliance with ISO Guide 35 & ISO/IEC: 17025.

Specification and method used Indian pharmacopoeia.
CRM/RM lot IPRS/M/028 is used for the comparison.

COMMUTABILITY: Not Applicable

INTENDED USE:

PILTRM/CRM is intended for use in product/material testing/calibration including R&D, Validation or Quality Control of Analytical Methods with specified quantity. This Material cannot be used as "Drug" or household.

INSTRUCTION FOR HANDLING & USE:

Allow the sealed container to equilibrate at room temperature before opening for use. Do not dry, use "On as is Basis". Once the container has been opened, Stability of content, value cannot be guaranteed. It is for immediate use. Read MSDS before use.



PRIYAN INTERNATIONAL LAB AND TECHNOLOGY

(Reference Material Product as per International Standard ISO: 17034)

2nd and 3rd Floor, C-247, Sector-10, Noida, Gautam Buddha Nagar, Uttar Pradesh-201301

Tel No. 0120-3684527, Mob.: 8882764797

Mail: priyanintlabetech@gmail.com, Web: www.priyaninternationalallabtech.com

CRM CERTIFICATE

Format No: PILT/QSP/055/00/FMT/02

VALIDITY:

Stated Validity is apply, when material stored under recommended conditions with proper handling. Any change in assigned value due to stability/retesting/review etc. or validity extension/revalidation/Updates, will be made available on our Website:

www.priyaninternationalallabtech.com

SAFETY INFORMATION:

Refer to the material safety data sheet.

Sandeep
20/06/2024

Approving Authority


MATERIAL SAFETY DATA SHEET (MSDS)

Company Information

Name of organization : PRIYAN INTERNATIONAL LAB AND TECHNOLOGY
Address : C-247, 2nd & 3rd Floor, Sector-10, Noida-201301
Ph. No. : 0120-3684527, +91-8882764797
Email : privanintlabbtech@gmail.com
Website : www.privaninternationalabbtech.com

Section-1. Product Identification and Composition

Product Name	: Methyl Cobalamin
Product No.	: M-009
Uses	: Laboratory chemicals, Reference Material
Chemical Formula	: $C_{63}H_{91}CON_{13}O_{14}P$
CAS No.	: 13422-55-4
Molecular Weight	: 1344.40 g/mol



Section-2. Hazards Identification

Potential Acute Health Effects:

Hazardous in case of skin contact (irritant), of eye contact (irritant), of ingestion, of inhalation. Severe over-exposure can result in death.

Potential Chronic Health Effects:

CARCINOGENIC EFFECTS: 3 (Not classifiable for human.) by IARC. **MUTAGENIC EFFECTS:** Mutagenic for mammalian somatic cells. Mutagenic for bacteria and/or yeast. **TERATOGENIC EFFECTS:** Not available. **DEVELOPMENTAL TOXICITY:** Not available. The substance may be toxic to heart, gastrointestinal tract, central nervous system (CNS). Repeated or prolonged exposure to the substance can produce target organs damage. Repeated exposure to a highly toxic material may produce general deterioration of health by an accumulation in one or many human organs.

Section-3. First Aid Measures

Ingestion:

Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.

Inhalation:

If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician

Skin Contact:

Wash off with soap and plenty of water. Consult a physician.

Eye Contact:

Flush eyes with water as a precaution.

Section-4. Fire and Explosion Data

Flammability:

May be combustible at high temperature.

Flammable Limits:

Not available.

Auto-Ignition Temperature:

Not available.

Extinguishing media

- ❖ Suitable extinguishing media - Use water spray, alcohol- resistant foam, dry chemical or carbon dioxide.
- ❖ Unsuitable extinguishing media - No data available.

Protection against fire:

Wear suitable protective equipment.

Hazardous combustion products:

No data available.

Section-5. Accidental Release Measures

Environmental precautions:

Prevent further leakage or spillage if safe to do so. Do not let product enter drains. Discharge into the environment must be avoided.

Personal precautions, protective equipment and emergency procedures:

Use personal protective equipment. Avoid dust formation. Avoid breathing vapors, mist or gas.

Ensure adequate ventilation. Avoid breathing dust.

Methods and materials for containment and cleaning up:

Pick up and arrange disposal without creating dust. Sweep up and shovel. Keep in suitable, closed containers for disposal.

Section-6. Handling and Storage

Hygiene measures:

Handle in accordance with good industrial hygiene and safety practice. Wash hands before breaks and at the end of workday.

Handling:

Avoid contact with skin and eyes. Avoid formation of dust and aerosols. Provide appropriate exhaust ventilation at places where dust is formed.

Storage:

Keep container tightly closed. Keep container in a cool, well-ventilated area.

Section-7. Exposure Controls/ Personal Protection

Respiratory Protection:

For nuisance exposures use type P95 (US) or type P1 (EU EN 143) particle respirator. For higher level protection use type OV/AG/P99 (US) or type ABEK-P2 (EU EN 143) respirator cartridges. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Hand protection:

Wear nitrile or other impervious gloves if skin contact is possible. When the material is dissolved or suspended in an organic solvent, wear gloves that provide protection against the solvent.

Skin Protection:

Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.

Eye Protection:

Safety glasses with side-shields conforming to EN166 Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

Body protection:

Complete suit protecting against chemicals, the type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

General Hygiene Consideration:

Prevent further leakage or spillage if safe to do so. Do not let product enter drains. Discharge into the environment must be avoided.

Section-8. Physical and Chemical Properties**8.1 Information on basic physical and chemical properties**

- | | |
|---|--|
| a) Physical state | : Crystalline powder |
| b) Color | : No data available |
| c) Odor | : No data available |
| d) Melting point/freezing point | : No data available |
| e) Initial boiling point and boiling range | : No data available |
| f) Flammability (solid, gas) | : No data available |
| g) Upper/lower flammability or explosive limits | : No data available |
| h) Flash point | : No data available |
| i) Auto ignition temperature | : No data available |
| j) Decomposition temperature | : No data available |
| k) pH | : No data available |
| l) Viscosity | : Viscosity, kinematic: No data available
: Viscosity, dynamic: No data available |
| m) Water solubility | : No data available |
| n) Partition coefficient n-octanol/water | : No data available |
| o) Vapor pressure | : No data available |
| p) Density | : No data available |
| Relative density | : No data available |
| q) Relative vapor density | |
| r) Particle characteristics | : No data available |
| s) Explosive properties | : No data available |
| t) Oxidizing properties | : none |

Section-9. Stability and Reactivity

Reactivity:

No data available

Chemical stability:

Stable under recommended storage conditions.

Possibility of hazardous reactions:

No data available

Conditions to avoid:

No data available

Incompatible materials:

Strong oxidizing agents

Section-10. Toxicological Information

Acute toxicity

LD50 Oral - Rat - > 5.000 mg/kg

Inhalation: No data available

Dermal: No data available

Skin corrosion/irritation

No data available

Serious eye damage/eye irritation

No data available

Respiratory or skin sensitization

No data available

Germ cell mutagenicity

No data available

Carcinogenicity

No data available

Reproductive toxicity

No data available

Specific target organ toxicity - single exposure

No data available

Specific target organ toxicity - repeated exposure

No data available

Aspiration hazard

No data available

Section-11. Ecological Information

Toxicity:

No data available

Persistence and degradability:

No data available

Bioaccumulative potential:

No data available

Mobility in soil:

No data available

Results of PBT and vPvB assessment:

This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

Section-12 Disposal Considerations

Waste treatment methods

Product

Offer surplus and non-recyclable solutions to a licensed disposal company. Dissolve or mix the material with a combustible solvent and burn in a chemical incinerator equipped with an afterburner and scrubber. Waste material must be disposed of in accordance with the Directive on waste 2008/98/EC as well as other national and local regulations. Leave chemicals in original containers. No mixing with other waste. Handle uncleaned containers like the product itself.

Contaminated packaging

Dispose of as unused product.

Section-13. Transport Information

UN number

ADR/RID: -

IMDG: -

IATA: -

UN proper shipping name

ADR/RID: Not dangerous goods

IMDG: Not dangerous goods

IATA: Not dangerous goods

Transport hazard class(es)

ADR/RID: -

IMDG: -

IATA: -

Packaging group

ADR/RID: -

IMDG: -

IATA: -

Environmental hazards

ADR/RID: no

IMDG Marine pollutant: no

IATA: no

Special precautions for user

Further information

Not classified as dangerous in the meaning of transport regulations.

Section-14. Regulatory Information

Safety, health and environmental regulations/legislation specific for the substance or mixture

This material safety data sheet complies with the requirements of Regulation (EC) No.1907/2006.

Chemical Safety Assessment

For this product a chemical safety assessment was not carried out.

Section-15. Other Information

Important Notice

Information applies only to this material of its intended use.

The PILT prepares the MSDS by using information available at the time from sources considerable, reliable, such as PILT approved summaries of product characteristics, RTECS and the MSDS of the suppliers, manufacturers or importers. The PILT does not independently verify the information. The accuracy of the information can't therefore be guaranteed, nor does it constitute any expression of opinion by the PILT concerning the Reference Material preparation. This information is accordingly not to be regarded as a representation or statement concerning the quality or safety of the Reference Material, the presence of any defect in it, or its fitness for any particular purpose except that of use as a IPRS by professional persons having technical skill and at their own discretion and risk. The downstream users have the responsibility to manage the risks arising from their usage of the PILT Reference Material and for use of any information provided in this MSDS. People working with any reference material should apply regional and national laws, good practices and state of the art precautions.